



Reewarm™

PTA Balloon
Dilatation Catheter


*Through
Anyway*





Composite Hydrophilic Coating

Improve the performance of passing capacity





Inner Hypo Tube with Metal Cutting

The length of spiral cutting is corresponding to the length of push pole, which ensure the excellent deliver ability and trackability to meet different clinical needs



Unique Ballon Material

Reduce the damage to blood vessel and the risk of dissection

Balloon Diameter (mm)	Balloon Length (mm)	Catheter Length (mm)	NP (atm)	RBP (atm)	Guide-wire (inch)	Catheter Type	Balloon Material	Pleat and Fold
2.0	20/40/60/ 80/100/ 120/150/ 220	90/150	6	14	0.018''	OTW	Nylon	3-Folding 
2.5		90/150	6	14				
3.0		90/150	6	14				
3.5		90/150	6	14				
4.0		90/130/150	6	14				6-Folding 
5.0		90/130	6	13				
6.0		90/130	6	13				
7.0		90/130	6	12				



Shanghai MicroPort Endovascular MedTech Co., Ltd.

Address: Building #1, 3399 Kangxin Rd., SIMZ Century Medicine Park, Shanghai 201318, P. R. China

Tel: +86 21 38139300

Fax: +86 21 33750026

Website: www.endovastec.com

All cited trademarks are the property of their respective owners in respective countries. These devices shall only be used for curing diseases on prescription by a licensed physician. Indications, contraindications, warnings, potential adverse events and directions for use can be found in Instructions for Use and product labeling supplied with each device. These devices are prohibited to promote, distribute and sell in any countries where the registration certificates or any other licenses are not obtained. This brochure is not intended for use in aforesaid countries. This brochure is only delivered to specific people for their reference, is not to be distributed widely as an advertisement. Any information may be changed without prior notice.

Not Available for Sale in the United States.

©2019 Shanghai MicroPort Endovascular MedTech Co., Ltd. All rights reserved. Version No. IB-RW-19-01